



Stainless Steel Bloodless Burdizzo Castrator Veterinary Instruments Simple Operation

Our Product Introduction

Basic Information

- Place of Origin: Jiangsu China(Mainland)
- Brand Name: CHUANGPU
- Certification: ISO
- Model Number: HL-Q2A Q2, Q2B ,Q2C
- Minimum Order Quantity: 1set
- Price: USD\$95.00/set
- Packaging Details: CARTON
- Payment Terms: Western Union, MoneyGram, T/T, L/C
- Supply Ability: 100SET/MONTH



Product Specification

- Features: Simple Operation Easy To Use Bloodless Castration
- Material: Stainless Steel
- Application: Horse, Cow, Goat, Pig
- Condition: New
- Highlight: **Stainless Steel Bloodless Burdizzo Castrator, Veterinary Instruments Bloodless Burdizzo Castrator**
,
Simple operation Bloodless Burdizzo Castrator

Product Description

Burdizzo castrator veterinary instruments , cow bull castrator emasculator, goat and sheep castrating tools

Calves must be castrated as young as practically possible, and pain control is required when castrating bulls older than six months of age. The requirement to use pain control in older calves was based on research demonstrating its effectiveness in feedlot bulls.

Sheep and goat breeders wait until at least two months to castrate because it gives the urinary tract more time to mature. This is especially important for kids that will be pets or weed eaters, as well as fiber wethers because those animals will hopefully live long lives.

A gelding is a castrated horse or other equine, such as a donkey or a mule. Castration, as well as the elimination of hormonally-driven behavior associated with a stallion, allows a male horse to be calmer and better-behaved, making the animal quieter, gentler and potentially more suitable as an everyday working animal.

Most pig producers castrated their pigs **between day four and 14 of age**. If a producer waits longer than 14 days, AVMA recommends that pigs be castrated at least **five days before weaning**.

Model	HL-Q2	HL-Q2A	HL-Q2B	HL-Q2C
Name	Bloodless Castrating(Burdizzo)			
Material	Stainless steel			
Animal	Cow	Horse, Deer, Mule	Pig / Buck	Goat
Weight	2.18kg	1kg	0.65kg	1.8kg



产品简介

Product introduction

创普机械 Chuang Pu

Product name: Bloodless Castrator

Product introduction: The bloodless castrator is a scientific veterinary surgical instrument. It is suitable for castration of livestock such as bulls, goat and pigs.

The device achieves the purpose of surgery by clamping the animal spermatic cord through the scrotum of the livestock, and does not need to incision on the scrotum of the livestock, so it is called "no blood castration".

Bloodless Castrator for cow



Bloodless castrator for goat



Bloodless castrator for pig




Bloodless castrator for horse, deer and Mule



Using bloodless castration method: livestock binding (veterinary terminology, the reliable and stable control of livestock), after the local anesthesia, the operator grasps the scrotum neck of livestock (such as cattle, sheep, etc.) by hand, and squeezes the testicles to the bottom of the scrotum, the spermatic cord is pushed to the outside of the neck of the scrotum, and is clamped on the inner skin of the spermatic cord with a long-handled sperm stringing clamp to prevent the spermatic cord from sliding and then open the jaw of the bloodless castration clamp and clamp it in the long handle. The ejaculation fixed clamp is fixed at 3-5CM above the fixed point, and the assistant slowly closes the handle. The clamp point should be at least 1 cm above the testicle.



After confirming that the spermatic cord has been clamped by the jaws, forcefully close the gripper handle and you will hear a crisp "click" sound, indicating that the spermatic cord has been chipped. After the clamp handle is closed, it should be kept for at least one minute, then loosen the jaws to ensure that the spermatic cord has been broken, loosen the forceps, and clamp the second time on the spermatic cord at 1.5-2.0CM below it to ensure the surgical effect.



The left picture is a schematic diagram of the surgical procedure



1. The gap between the jaws is tight, the success rate of the operation is high.



Hailian Packaging Equipment Co.,Ltd

☎ 0086 13301569240

✉ sales2@dairymachinery.cc

🌐 mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China