

## Sheep Bloodless Castrator Stainless Steel Cultivation Farm Using Tool

Our Product Introduction

for more products please visit us on [mobilemilkingmachine.com](http://mobilemilkingmachine.com)

### Basic Information

- Place of Origin: Jiangsu
- Brand Name: chuangu
- Certification: CE ISO FDA SGS
- Model Number: HL-Q2C
- Minimum Order Quantity: 1
- Price: USD\$30/set
- Packaging Details: Carton
- Delivery Time: 5-8days
- Payment Terms: T/T, Western Union,Paypal
- Supply Ability: 100set/month




### Product Specification

- Material: Stainless Steel
- Package: Carton
- Usage: Sheep Castrator
- Model: HL-Q2C
- MOQ: 1 PCS
- Delivery Time: 2-3days
- Highlight: cattle hoof trimming discs, cattle nose ring



### More Images



## Product Description

### Sheep Bloodless Castrator Stainless Steel Cultivation Farm Using Tool

#### Product Introduction:

Product name: Sheep Bloodless Castrator

Product number: HL-Q2C

Product introduction: The reason we then called "no blood castration" is the bloodless castrator is a scientific veterinary surgical instrument. The device achieves the purpose of surgery by clamping the animal spermatic cord through the scrotum of the livestock, and does not need to incision on the scrotum of the livestock. It is suitable for castration of livestock such as bulls, sheep and pigs. Easy to cooperation and Any question you can contact us any time.

#### Bloodless castrator for goat



Model	HL-Q2C
Material	Stainless Steel
Package	Carton
Usage	Sheep Castrator
Animal	Sheep or Goat
Weight	1.14KG
Country	China
Delivery time	3-5 days

#### Using bloodless castration method:

livestock binding (veterinary terminology, the reliable and stable control of livestock), after the local anesthesia, the operator grasps the scrotum neck of livestock (such as cattle, sheep, etc.) by hand, and squeezes the testicles to the bottom of the scrotum, the spermatic cord is pushed to the outside of the neck of the scrotum, and is clamped on the inner skin of the spermatic cord with a long-handed sperm stringing clamp to prevent the spermatic cord from sliding and then open the jaw of the bloodless castration clamp and clamp it in the long handle. The ejaculation fixed clamp is fixed at 3-5CM above the fixed point, and the assistant slowly closes the handle. The clamp point should be at least 1 cm above the testicle. After confirming that the spermatic cord has been clamped by the jaws, forcefully close the gripper handle and you will hear a crisp click" sound, indicating that the spermatic cord has been clipped. After the clamp handle is closed, it should be kept for at least one minute, then loosen the jaws to ensure that the spermatic cord has been broken, loosen the forceps, and clamp the second time on the spermatic cord at 1.5-2.0CM below it to ensure the surgical effect.



**Hailian Packaging Equipment Co.,Ltd**

☎ 0086 13301569240

✉ sales2@dairymachinery.cc

🌐 mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China