



Durable Teeth Cutting Pliers For Rabbit , Stainless Steel Pig Teeth Clipping

Basic Information

Place of Origin: china
Brand Name: chuangpu
Certification: ISO
Model Number: HL-Q17
Minimum Order Quantity: 1 unit

• Price: USD\$4.00/unit

Packaging Details: Carton

• Delivery Time: 7-10 working days

Payment Terms: T/T, Western Union, MoneyGram,Paypal

• Supply Ability: 1000 units per month



Product Specification

Name: Teeth Plier
Material: Stainless Steel
Animal: Pig And Rabbit
Function: Cutting Teeth
Model: HL-Q17

• Highlight: cattle hoof trimming discs, cattle nose ring



More Images









Product Description

Pig tooth cutting plier, teeth cutting pliers for rabbit, stainless steel pig teeth clipping

The **teeth plier** is stainless steel **pig teeth** clipping. This **teeth** clipping is only suitable for **cutting** small **teeth** such as piglets and rabbits.

Soon after they are born, the **teeth** of piglets are often clipped. The purpose of **teeth**-clipping is to reduce injuries caused to each other and to their mother as piglets nurse.

Model	HL-Q17
Name	Teeth plier
Material	Stainless Steel
Animal	Pig / Piglet, rabbit
Function	cutting / clipping teeth



创普机械 Chuang Pu



Why use a forceps







- 1. When piglets compete for nipples, they use teeth to fight each other and bite their cheeks. Once infected by bacteria, they will affect their milk consumption, and some even cause blinking and loss of value.
- 2. When the piglets compete for the nipple, the sharp canine teeth will bite the sow's nipple or breast, causing breast inflammation and sow pain, refuse to breastfeed, seriously affecting the growth and development of the piglet.
- 3. Sows often lie down due to pain, stepping on piglets
- 4.Pigs tend to bite during the period of conservation and fattening, which seriously affects breeding

创普机械 Chuang Pu

Instructions



- 1.Pinch the upper and lower jaws of the piglets with your thumb and forefinger, forcing the piglets to open their mouths to reveal their teeth, and then use the diagonal pliers of the forceps to cut the left and right teeth. Be careful not to injure gums and tongue when cutting teeth. The cut tooth granules do not allow the piglets to swallow.
- 2.PEvery tooth cut of a piglet should be disinfected by the tooth cutting tool to avoid cross-infection of the bacteria; if the piglet is weak after birth, the tooth will not be cut temporarily to strengthen its ability to compete for the nipple, and then cut the tooth after 3–4 days.
- 3.After the piglet is born, 8 sharp needle-shaped teeth are cut off, which can reduce the damage to the sow's nipple. When the battle occurs, the damage to the litter can also be reduced.





Hailian Packaging Equipment Co.,Ltd

0086 13301569240

sales2@dairymachinery.cc

mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China