

Jiangsu

HL-Q3

Carton

2-4 days

Chuangpu

CE ISO FDA SGS

USD\$11-12/pcs

500PCS/Month

# Metal Veterinarian Tools Bull - Holder Bull Nose Tong Cattle Head Holder

# **Basic Information**

- Place of Origin:
- Brand Name:
- Certification:
- Model Number:
- Minimum Order Quantity: 1pcs
- Price:
- Packaging Details:
- Delivery Time:
- Payment Terms: T/T, Western Union, Paypal
- Supply Ability:



রাগ

### **Product Specification**

- Material:
- Package:
- Usage:
- Length:
- Weight:
- Delivery Time:
- Highlight:

Carbon Steel
Carton
Bull Holder
350mm
0.77KG
2-4 Days

cattle hoof trimming discs, cattle nose ring



### More Images



#### Metal Veterinarian Tools Bull - Holder Bull Nose Tong Cattle Head Holder

## Poduct Information

### Product name: Cattle nose pliers

Product number: HL-Q3,Q4,Q4A,Q4C

Product introduction: The cattle nose plier is also called the cow nose clip. It is a kind of bull head sealer. When you use it, you only need to put the nose clamp on both sides of the nasal septum to transfer the attention of the cow, so that the cow can be injected intravenously to restrain the cow. There is so many style you cna choose for cow holder, any question you can contact me any time.

Model	HL-Q3
Material	Stainless Steel
Package	Carton
Usage	Bull Holder
Animal	Cow
Weight	0.77KG
Country	China
Delivery time	3-5 days





Product number: HL-Q3 Product material: Stainless steel Product specification: length35cm\ maximum opening16cm \ weight0.77kg Product Usage: Fixed cattle nose

# How to Use the Nose pliers?

1.prepare the material, take thick thumb rope, 1.5 meters long, a pair of cattle nose pliers.

2.Put a slip knot on the end of the rope, and fix it on the neck of the pliers so as not to affect the opening of the pliers.

3. When clamping the cattle nose with pliers, tighten the small rope, and then put a slip knot to fix the knot on the end of the two handles of the pliers.

4. Tighten the rope and tie the rope to the pile.



