

Pill Propellant Gun Dairy Machinery Appliance Dairy Farm Animal Feeding Pill Tool

Our Product Introduction

for more products please visit us on mobilemilkingmachine.com

Basic Information

- Place of Origin: jiangsu
- Brand Name: chuangu
- Certification: ISO CE FDA SGS
- Model Number: HL-Q11
- Minimum Order Quantity: 100pcs
- Price: USD\$20-25/pcs
- Packaging Details: carton
- Delivery Time: 2-4days
- Payment Terms: T/T, Western Union,Paypal
- Supply Ability: 2000pcs/month



Product Specification

- Material: Stainless Steel
- Length: 650mm/700mm
- Package: Carton
- Weight: 0.5kg/1.15kg
- Usage: Pill Feeding
- Style: 2 Styles
- Highlight: cattle hoof trimming discs, cattle nose ring



Product Description

Product introduction

The pill propellant gun is made of stainless steel material and is used to administer the medicine when the animal is sick. It is a simple and convenient instrument for cattle and sheep. In operation, put the medicine to be fed to the animal into the medicine bucket and then press the awl firmly. After the completion of the filling, the injection gun should be cleaned in time for the next use.

Product information

Model: HL-Q11 & HL-Q11A

Item: pill propellant gun

Usage: feeding animal

Package: Carton

Length: 650/700mm

Model	HL-Q11	HL-Q11A
Material	stainless steel	stainless steel
Length	650mm	700mm
Package	carton	carton
Weight	0.5kg	1.15kg
Usage	animal feeding	animal feeding

产品参数

Product parameters

创普机械 Chuang Pu

款式一

Style 1



款式二

Style 2



Object that can be placed

创普机械
Chuang Pu



Magnetic Trap



Veterinary drug

Instructions

- 1.Put the drug or magnetic cage that needs to be fed into the head of the vat.
- 2.Open the animal' s mouth with an opener and then insert the medicine gun into the animal's mouth.
- 3.After confirming the successful investment, the medicine gun can be taken out.

Product details



Hailian Packaging Equipment Co.,Ltd

☎ 0086 13301569240

✉ sales2@dairymachinery.cc

🌐 mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China