

## Livestock Pig Cattle Sheep Thawing Cup Aluminum Alloy ABS For Optimal Thawing

Our Product Introduction

for more products please visit us on [mobilemilkingmachine.com](http://mobilemilkingmachine.com)

### Basic Information

- Place of Origin: Jiangsu, China(Mainland)
- Brand Name: Chuangpu
- Model Number: HL-FS-TTC
- Minimum Order Quantity: 1 PCS
- Price: Negotiable
- Packaging Details: Carton
- Delivery Time: 7-10 working days
- Payment Terms: T/T, Western Union,
- Supply Ability: 5000PCS /Month



### Product Specification

- Name: Frozen Semen Thawing Thermostatic Cup
- Material: Aluminum Alloy+ABS
- Volume: 600ml
- Applicable Objects: Cattle, Sheep, Pigs,
- Weight: 1.4KG
- Highlight: **Optimal Thawing Sheep Thawing Cup, Aluminum Alloy Sheep Thawing Cup, Livestock Cattle Sheep Thawing Cup**



### More Images



## Product Description

**Corporate HQ:** [CHUANGPU ANIMAL HUSBANDRY TECHNOLOGY\(SUZHOU\) CO.,LTD](#)

**Export Dept:** [HAILIAN PACKAGING EQUIPMENT CO., LTD](#)

**Branch:** [Sichuan Gene Haven Animal Husbandry Technology Co.Ltd](#)

EMAIL: [nancyzd68@gmail.com](mailto:nancyzd68@gmail.com)  
[sales2@dairymachinery.cc](mailto:sales2@dairymachinery.cc)

MOBILE: +86 133 0156 9240(Wechat/WhatsApp)  
 OFFICE: +86 512 5890 5953

### Semen Thawing Cup, Ideal Frozen Sperm Thawing Cup for Livestock Pig Cattle Sheep

The frozen semen thawing cup is primarily used to thaw frozen semen for artificial insemination

The frozen semen thawing cup is a specialized device for thawing frozen semen. It provides precise temperature control, ensuring that the semen is not damaged during the thawing process. When using, simply place the frozen semen into the thawing cup and wait for a while to complete the thawing process. This kind of device usually has a portable design, suitable for performing artificial insemination operations in different locations. In addition, some high-end thawing cups are equipped with timers and other smart functions to better control thawing time and temperature.

Item Type: Frozen Thawing Cup  
 Material: Aluminum Alloy+ABS  
 Weight: Approx. 1358g / 47.9oz  
 Product Size: Approx. 21x10.5cm / 8.3x4.1in  
 Applicable Objects: pigs, cattle, sheep, dogs  
 Voltage Specifications: 100-240V US Plug  
 Capacity: 600ML  
 Battery (Built-In Battery): 1.1V lithium battery, 13200mAh capacity  
 Input Voltage: 12-14V  
 Input Current: 5-7A  
 Charging Time: 3-5 hours  
 Refreshing Frequency: 2 (times/second)  
 Use Environment: 0-40  
 Buffer Temperature Difference: 0.1-0.3 degrees Celsius  
 Standby at 20°C Constant Temperature: 4-5 hours  
 Standby at 0°C Constant Temperature: 2 hours  
 Standby at -20°C Constant Temperature: 40 minutes  
 Introduction: is highly sensitive to temperature. D

#### Feature:

1. The bottom of the cup has a non-slip mat design. Multi-temperature adjustable, used for thawing frozen and constant temperature of sperm.
2. A variety of power supply methods can be used with vehicle power cords and household power supplies.
3. 600ML aluminum alloy inner, covering heating, rapid and accurate heating.
4. Touch-type numerical control operation, multi-temperature adjustment and convenient operation. The specified time can be set to achieve constant temperature heating within the specified time.
5. The built-in storage battery can run independently to keep the temperature in the cup constant at the set temperature and keep the semen in the best condition. Dilute it at the minimum temperature difference, which helps to extend the storage time of the sperm.



Switch on and off: after the power is turned on or the switch is turned on, the key background light lights up, long press the on key to start the machine, and long press the on key to shut it off under the starting state.

Temperature setting: press the setting button, hear the contact sound, the temperature value flashes, adjust the temperature by adding or subtracting on the right side, and automatically confirm 10 seconds later.

Timer setting: press the setting button twice, hear two electric shocks, the timer flashes, adjust the temperature by addition or subtraction, automatically confirm 10 seconds later, press the open key to start the timer, and the alarm will sound when the time is up.



0086 13301569240



[sales2@dairymachinery.cc](mailto:sales2@dairymachinery.cc)



[mobilemilkingmachine.com](http://mobilemilkingmachine.com)

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China