

Cow Ruminal Puncture,Rumen Vent Needle,Stainless Steel Puncture Needle For Cattle Trocar And Canula Treatment ,

Our Product Introduction

for more products please visit us on mobilemilkingmachine.com

Basic Information

- Place of Origin: Jiangsu China(Mainland)
- Brand Name: CHUANGPU
- Certification: ISO
- Model Number: HL-Q8A
- Minimum Order Quantity: 1set
- Price: USD\$4.00/SET
- Packaging Details: CARTON
- Payment Terms: Western Union, MoneyGram, T/T, L/C
- Supply Ability: 1000SET/MONTH



Product Specification

- Material: Food Grade Stainless Steel
- Item: Rumen Vent Needle
- Farm Animals Type: Cattle, Cow
- Length: 200mm
- Usage: Cow Ruminal Puncture
- Highlight: cattle hoof trimming discs, cattle nose ring



More Images



Product Description

Cow Ruminal Puncture,Rumen Vent Needle,Stainless Steel Puncture Needle For Cattle Trocar And Canula Treatment ,

The rumen is the first chamber of the ruminant stomach. It is the largest chamber and has regular contractions to move food around for digestion, eliminate gases through eructation and send food particles back to the mouth for remastication. The rumen breaks down food particles through mechanical digestion and fermentation with the help of symbiotic microbes. Volatile fatty acids are the main product of ruminant digestion.

Product Model Number: HL-Q8A

Product Material: Food grade stainless steel

Product Specification: Length is 20mm, Net Weight is 0.08Kg.

Product Description: Using the general method of urgency, rumen puncture.

Exhaust method: the use of trocar needle in the rumen parts puncture deflation,

Deflation should be the first fast after the slow, intermittent deflation.

In the process of deflation should pay attention to not excessive, or risk caused by collapse.

Farm Animals Type	Cattle/Cow/Horse/Camel etc
Item	Cow Ruminal Puncture

Our Product Introduction

Bleed Needle

HOT

产品简介 Product Introduction

创普机械 Chuang Pu

创普机械
Chuang Pu



Product name: Bleed Needle

Product number: HL-Q8A

Product introduction: It is used to vent the rumen of cattle, when the rumen bulging of cattle, sheep, camels, etc., the general first aid method, rumen puncture method (exhaust method)

产品参数 Product parameters

创普机械 Chuang Pu

High quality stainless steel



Weight: 0.08KG



创普机械
Chuang Pu

Needle sleeve

Bleed Needle

Instructions for use

Using a trocar to puncture and deflate the rumen, the deflation should be fast and slow, with intermittent deflation. In the process of deflation, care should be taken not to overdo it, otherwise there is a danger of causing collapse.



创普机械
Chuang Pu



1. Made of high quality stainless steel with sharp needles.



2. Equipped with a needle sleeve to avoid injury and convenient storage.



创普机械
Chuang Pu

3. High quality plastic grip for easy access.



E-mail: sales2@dairymachinery.cc
nancyzxd68@gmail.com

Skype: nancy.zhang62
Mobile: +86 133 0156 9240 (Wechat/WhatsApp)



Hailian Packaging Equipment Co.,Ltd



0086 13301569240



sales2@dairymachinery.cc



mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China