



Ruminal Magnet, Strong Rumen Magnet, Cow Magnet, strong magnetic cage rumen iron absorption

Our Product Introduction

for more products please visit us on mobilemilkingmachine.com

Basic Information

- Place of Origin: Jiangsu, China(Mainland)
- Brand Name: Chuangpu
- Certification: ISO
- Model Number: HL-Q12
- Minimum Order Quantity: 1pcs
- Price: USD\$18.00/PCS-USD\$18.50/PCS
- Packaging Details: CARTON
- Payment Terms: T/T, L/C, Western Union, MoneyGram
- Supply Ability: 1000PCS PER MONTH



Product Specification

- Features: Simple Operation Easy To Use
- Brand Name: Chuangpu
- Model Number: HL-Q12
- Animal Type Application: Cow/Cattle
- Place Of Original: Jiangsu China(Mainland)
- Is_customized: Yes
- Condition: New
- Item: Take Iron Wire Or Iron Nail Form Cow Stomach Machine
- Material: Stainless Steel
- Highlight: veterinary medical instruments, cattle hoof trimming discs



More Images



Product Description

Ruminal Magnet, Strong Rumen Magnet, Cow Magnet, strong magnetic cage rumen iron absorption

A rumen magnet is like a very large pill that that animal is forced to swallow with a pill gun. The magnet stays in the cow's rumen for the life of the animal and traps any metal the animal may swallow.

A rumen magnet attracts metal objects to help them pass through your cow's digestive system without causing internal injury.

Rumen magnets help cattle farmers prevent hardware disease in the herd. When cows eat small pieces of metal such as screws and wire, it can lead to numerous digestive issues and other health problems.



产品简介

Product Introduction

创普机械 Chuang Pu




Product name: Magnet Absorbers

Product number: HL-Q12

Product introduction: It is a tool for removing nails, iron wire and other iron scraps from the stomach of cattle. It can effectively prevent diseases such as traumatic net gastritis, pericarditis and pleurisy in cattle, thereby reducing the mortality of cattle.

Feeder



It consists of three tubular rods (stainless steel) that are threadedly connected to each other, and a flexible sleeve at the end of the ram is convenient for hand-held.

Holder



A foldable retaining ring is attached to each end of a main frame. The retaining ring is connected with a belt for tethering

创普机械
Chuang Pu

Absorbers

A foldable retaining ring is attached to each end of a main frame. The retaining ring is connected with a belt for tethering. It can absorb objects with a gravity of not less than 5 kg on the surface of 4 cm square of the contact surface, and the magnetic force can be maintained for 3 years.




Instructions

Open the cow's mouth, put the holder into the mouth of the cow (the square hole is external), and fix the fixing belt to the back side of the two corners of the cow. Then, the magnetic head with the lead-out rope of the feeder is inserted into the beef stomach through the square hole of the fixer, the feeder is held by one hand, and the lead rope is grasped by one hand, and the take-up rope is gently tensioned together with the feeder when taking out. Finally remove the holder.



1. The rear end of the magnet block shell has an elliptical shape to prevent the esophagus from .



2. The push rod is divided into three sections, which are easy to assemble.



3. Widened plastic handle for easy grip.



Hailian Packaging Equipment Co.,Ltd

☎ 0086 13301569240

✉ sales2@dairymachinery.cc

🌐 mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China