



# Animal Care Cow Drencher Pneumatic Cylinder Stomach Pump System With Nose Clip

### **Basic Information**

• Place of Origin: Jiangsu China(Mainland)

Brand Name: CHUANGPU

Certification: ISOModel Number: HL-Q1Minimum Order Quantity: 1set

• Price: USD\$300.00/set

Packaging Details: CARTON

• Payment Terms: Western Union, MoneyGram, T/T, L/C

• Supply Ability: 100SET/MONTH



# **Product Specification**

Material: Stainless Steel, PVC

Application: Cattle

• Item: Cow Rehydration Device Goat Fluid Infusion

Cattle Fluid Replacement

Form: Cattle Nose Clip Pneumatic Cylinder Hose

Sealing Ring Spring Hose

• Fluid Input Fluid Input 2-3min/50KG

Volume:

• Usage: Horse

Highlight: cattle hoof trimming discs, cattle nose ring



## More Images





#### Product Model Number: HL-Q1

### Animal care cow drencher, cattle stomach pump system, dairy cattle drenching with nose clip

The cow drencher is a unique pumping system specifically designed for high volume drenching pre or post fresh and sick cows and cattle.

It was designed for the sole purpose of making cow and cattle drenching easier, lowering the risk, taking less time to get the job done but

yet making it very affordable.

Product Material: Stainless Steel, PVC

#### Description:

Dairy rehydration device is composed of drug delivery device, transparent catheter and gastric tube three parts.

Wherein the drug pumping device through the transparent conduit tube link, trans-esophageal gastric tube is inserted into the rumen.

And then drug pumping device pumping liquid medicine delivery to the stomach tube, into the rumen, to complete the dehydration of cattle water.

#### Feature

Cattle Nose Clip: one hand operation, opening diameter of 70mm, for all types of cattle.

Pneumatic Cylinder: with a check valve to prevent liquid backflow.


Gas cylinders, handles, spring tube are made of 304 stainless steel and is equipped with non-slip rubber handle, cylinder pressure resistance does not exceed 1kg;

Hose: the use of food-grade transparent milk tube extrusion molding without seams, no leakage;

Sealing Ring: the use of nitrile rubber non-toxic life lasting.

Spring Hose: the top Seiko polished, no burr, chamfering the mouth of the cattle, the esophagus without any damage.

Name	Postpartum Rehydration Device
Model	HL-Q1
Material	Stainless steel
Consist of	Inflator, Soft Tube, Spring Hose, Bull-holder
Weight	5kgs
Function	administer fluids to a dehydrated cow







Product name: Cattle Fluid Feeder

Product number: HL-Q1

Product introduction: Cattle Fluid Feeder is a three-part instrument consisting of a drug pumping device, a transparent catheter and a stomach tube, wherein the drug pumping device is connected to the stomach tube through a transparent catheter, the stomach tube is inserted into the rumen of the cow through the esophagus, and then the medicine pumping device sucks. The liquid medicine is delivered to the stomach tube and enters the rumen to complete the hydration of dehydrated cattle.

使用说明 Directions for use

创普机械 Chuang Pu



# Prepare the filling solution



According to the appropriate proportion, flush the liquid with warm water and pour it into the bucket.



# Prepare the cows waiting to be fed



Keep the bull's head low and reach forward so that the oral esophagus is basically in a straight line through the rumen. Attach the tube conne nose clamp oral cavity.



Attach the esophageal guide tube connected with the bull nose clamp carefully to the oral cavity.



Using bull nose clamp to grip the cows' nose.



# Insert esophageal irrigation tube correctly Chuang Pu



Carefully insert the esophageal irrigation tube. If you feel that there is even resistance in the mouth of the cow, stop for a while and allow the cow to do his own swallowing and swallow the tube. Don't push it.



# Inject the liquid into the rumen



Put the pumping device into the bucket, press the pumping device by hand to pump the liquid medicine to the stomach tube, enter the rumen, and complete the hydration of dehydrated cattle.



### Hailian Packaging Equipment Co.,Ltd



0086 13301569240



sales2@dairymachinery.cc



@ mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China