

Castration Tool Castration Bander for Tail Docking Fits Bull,Calves, Goats, Sheep, Easy Application

Basic Information

• Place of Origin:

Model Number:

- Brand Name:
 - CHUANGPU HL-Q27B/HL-
- Minimum Order Quantity:
- Price:
- Packaging Details:
- Delivery Time:
- Payment Terms:
- Supply Ability:
- HL-Q27B/HL-Q27C 1 Bag Negotiable Carton 7-10 working days T/T, Western Union,

Jiangsu China(Mainland)

50000BAG /Month

Castration



রাগ

Product Specification

- Features:
- Material:
- Application:
- Condition:
- Highlight:

Special Pure Latex 、Silicone
Horse, Cow, Goat, Pig, Calf
New
Elastrator Bloodless Castrator,

Simple Operation Easy To Use Bloodless

Simple operation Bloodless Castrator, Simple operation bull castrator



More Images



Corporate HQ: CHUANGPU ANIMAL HUSBANDRY TECHNOLOGY(SUZHOU) CO., LTD

Export Dept: HAILIAN PACKAGING EQUIPMENT CO., LTD Branch: Sichuan Gene Haven Animal Husbandry Technology Co.Ltd

EMAIL: nancyzxd68@gmail.com sales2@dairymachinery.cc

MOBILE: +86 133 0156 9240(Wechat/WhatsApp) OFFICE: +86 512 5890 5953

Castration Tool Castration Bander for Tail Docking Fits Bull, Calves, Goats, Sheep, Easy Application

The tail ring is usually made of Special Pure Latex ,Silicone and is used by placing it around the base of the livestock's tail to hinder blood circulation, causing the tail to gradually shrink and eventually fall off This method is simple and effective, with no bleeding, fast healing, and good results.

There are two different castration rings for different animals.

Model:HL-Q27B Material:special pure latex Usage:Castrating or Cutting Tail Animal: lambs, piglets or other livestock Castration rings, Made of special pure latex material, non rubber, sturdy and strong toughness, durable and long service life

Model:HL-Q27C Material:Silicone Usage:Castrating or Cutting Tail Animal: calf Castration rings, Made of silicone material, sturdy and strong toughness, durable and long service life



Note: Pure manual measurement, if there is any error, please refer to the actual product

This type of tail docking and castration ring can be used for tail docking and animal castration in combination with tail docking forceps. It has the characteristics of simple postoperative treatment, no bleeding and no side effects.

Patien Packaging Equipment Co.,Ltd



sales2@dairymachinery.cc

@ mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China