



# 10mm Coarseness Bull Nose Ring With Durable Carbon Steel And Chrome Plating

# **Basic Information**

Place of Origin: China
Brand Name: Chuangpu
Model Number: HL-Q5D
Minimum Order 1 PCS

Quantity:

• Price: Negotiable

• Delivery Time: 3-5 working days

Carton

• Payment Terms: T/T

· Packaging Details:

• Supply Ability: 50000PCS /Month



# **Product Specification**

• Size: ID80mm; OD 100mm

Shape: Round

Auxiliary Tools: Bovine Nose Piercing Forceps Screwdriver

Material: Carbon Steel

Ring?:

Style: Reusable Bull Nose Ring

Weight: 0.17kgManufacturer: Chuangpu

Highlight: Durable Carbon Steel Bull Nose Ring,

Chrome Plating Bull Nose Ring,

10mm Bull Nose Ring



# More Images





# **Product Description**

Corporate HQ: CHUANGPU ANIMAL HUSBANDRY TECHNOLOGY(SUZHOU) CO.,LTD

Export Dept: ZHANGJIAGANG FANRS MACHINERY CO.,LTD. Branch: Sichuan Gene Haven Animal Husbandry Technology Co,Ltd

EMAIL: nancyzxd68@gmail.com sales2@dairymachinery.cc

MOBILE: +86 133 0156 9240(Wechat/WhatsApp)

OFFICE: +86 512 5890 5953

## **Product Description:**

The sliver color of this bull nose ring gives it a classic and timeless look that is sure to impress. It's perfect for those who are looking for a simple and elegant design that will look great on any animal. The round shape of this decorated nose ring gives it a unique and eyecatching appearance that will make your animal stand out from the crowd.

Included with this bull nose ring are the necessary auxiliary tools, including bovine nose piercing forceps and a screwdriver. This makes it easy to install and remove the nose ring without causing any unnecessary discomfort to your animal. The weight of this bull nose ring is 0.17kg, making it lightweight and easy to wear for your animal.

The Chuangpu Reusable Bull Nose Ring is a great investment for anyone looking for a high-quality animal nose ring that is both stylish and practical. Whether you're using it for shows, breeding, or just as a round on nose jewellery piece, this bull nose ring is sure to impress.

#### Features:

Product Name: Bull-Nose-Ring

Coarseness: 10mm Size: ID 70mm; OD 90mm

Shape: Round

Style: Reusable Bull Nose Ring Suitable For: Bull Nose Wear

This Animal Nose Ring is perfect for bull nose wear. The Bull Nose Ring is round and has a coarseness of 10mm, with a size of ID 80mm

and OD 100mm. The reusable style makes it practical and convenient for animal nose piercing.

## **Technical Parameters:**

Auxiliary Tools	Bovine Nose Piercing Forceps,Screwdriver
Suitable For	Bull Nose Wear
Durability	High
Shape	Round
Surface Treatment	Chrome Plating
Manufacturer	Chuangpu
Material	Carbon Steel
Size	ID 80mm; OD 100mm
Coarseness	10mm
Weight	0.17kg

# **Applications:**

Model Number: HL-Q5D Place of Origin: China Size: ID 80mm; OD 100mm

Auxiliary Tools: Bovine Nose Piercing Forceps, Screwdriver

Shape: Round Durability: High

Suitable For: Bull Nose Wear

Name: Chrome Cattle Nose Ring Product model: HL-Q5C, HL-Q5D

## Features:

1. Using chrome plating process, the surface is smooth and no burrs.

2. Self-locking design, convenient and practical.



Reminder: The above parameters are all measured manually, if there is any error, please refer to the actual product

## INSTALLATION METHOD

Step 1

Prepare the bull nose ring





Loosen the set screw with a screwdriver, then break the nose ring to the right size and prepare to put the ring on the cow.



Putting a bull's nose ring on a cow

Before putting on the nose ring, a person should hold the cow's neck and lift it upward. This way, the cow will not have the strength to resist. Then feel the thin area between the two nostrils of the cow's nose with your hand. Like the cartilage of a human ear. Put the bull nose ring through the bull nose and close it. Tighten the screw with the product, and then break the head of the screw with a little force.

### USING THE MATTERS NEEDING ATTENTION

- 1. After nose piercing, iodine tincture and antiseptic ointment can be applied to the wound (such as fish stone grease ointment).
- 2. The puncture point should not be too far in front, otherwise it will be easily broken; Too far back is not a good idea, because the nerves in the back are not sensitive and the nose ring is more difficult to penetrate.
- 3. After the nose ring is pierced, let the cattle rest for 1~2 days, and then often lead the cattle to a place with green grass to let them feed. Let the cattle know that the benefit of the nose ring is to have grass to eat, and gradually let the cattle do some farm work to train them to be docile.
- 4. Nose piercing should be performed at the first year of life (usually between 10 and 12 months of age).



Hailian Packaging Equipment Co.,Ltd







