


0.25ml Frozen Semen Tube 0.5ml Frozen Semen Straw Disposable PVC Clean Farming Accessory For Livestock

Our Product Introduction

for more products please visit us on mobilemilkingmachine.com


Basic Information

- Place of Origin: Imported from Belarus with original packaging
- Brand Name: Chuangpu
- Model Number: HL-FS-TT
- Minimum Order Quantity: 1 Bag
- Price: Negotiable
- Packaging Details: Carton
- Delivery Time: 10-15 working days
- Payment Terms: T/T, Western Union,
- Supply Ability: 50000BAG /Month




Product Specification

- Specification: 0.25ml And 0.5ml
- Length: 133mm
- Material: PVC
- 0.25ml OD: 1.92 - 1.96mm
- 0.25ml ID: 1.6 - 1.64mm
- 0.5ml OD: 2.87 - 2.91mm
- 0.5ml ID: 2.5 - 2.54mm
- Highlight: **PVC Frozen Semen Straw, Disposable Frozen Semen Straw, Farming Frozen Semen Tube**



More Images



Product Description

Corporate HQ: [CHUANGPU ANIMAL HUSBANDRY TECHNOLOGY\(SUZHOU\) CO.,LTD](#)

Export Dept: [HAILIAN PACKAGING EQUIPMENT CO., LTD](#)

Branch: [Sichuan Gene Haven Animal Husbandry Technology Co.Ltd](#)

EMAIL: nancyzxd68@gmail.com
sales2@dairymachinery.cc

MOBILE: +86 133 0156 9240(Wechat/WhatsApp)
OFFICE: +86 512 5890 5953

0.25ml Frozen Semen Tube 0.5ml Frozen Semen Straw Disposable PVC Clean Farming Accessory for Livestock

Freezing fine tubes are a standard tool widely used for freezing storage of semen from bulls, stallions, sires, and other species. The capacity of frozen fine tubes is typically 0.25ml or 0.5ml, making them ideal for artificial insemination.

These thin tubes are made from high-quality materials, ensuring that the quality of semen is not affected during freezing and thawing processes.

- **Premium PVC Material:** Made from premium PVC material, ensuring stability and reliability of the product.
- **High Capacity Design:** Each tube has a capacity of 0.25ml or 0.5ml, meeting the needs of animal semen freezing and providing ample capacity and convenience.
- **Disposable Use:** Designed for disposable use to avoid cross , simplify process, and ensure the safety and hygiene of the semen freezing process.
- **Standard Size:** Conforms to the common animal semen freezing plug on the market, making it convenient for replacement and application.
- **Complete Packaging:** Each pack contains 1000 tubes or 500 tubes, convenient for storage to avoid material loss and , ensuring the integrity and reliability of the tubes.

Specification:

Model: HL - FS - TT 0.25ml

Item Type: Frozen Semen Tube

Material: PVC

Application: Freeze semen

Capacity: 0.25ml

Length:133mm

OD: 1.92~1.96mm

ID: 1.6~1.64mm

Quantity: 1000PCS

Model: HL - FS - TT 0.5ml

Item Type: Frozen Semen Tube

Material: PVC

Application: Freeze semen

Capacity: 0.5ml

Length:133mm

OD: 2.87~2.91mm

ID: 2.5~2.54mm

Quantity: 500PCS

Why need use Frozen Semen Tube ?

The use of frozen fine tubes improves the utilization rate of excellent male livestock, accelerates breed improvement, enhances production performance, and is not limited by geography or time. It can significantly reduce the number of male livestock and save expenses.

Various colors are random and cannot be specified.



PRODUCT PARAMETERS

WE ALWAYS ADHERE TO RIGOROUS MANAGEMENT, THE FOUNDATION OF METICULOUS MANUFACTURING,
PRODUCE PRODUCTS THAT MEET THE MARKET'S RULES AND CUSTOMER NEEDS

Product name: Disposable veterinary frozen sperm tubes

Model: HL-FS-TT 0.25ml

Specifications: 0.25ml

Origin: Belarus

Raw Materials: Germany imported PVC

Packing: 1000 pieces/pack



Product name: Disposable veterinary frozen sperm tubes

Model: HL-FS-TT 0.5ml

Specifications: 0.5ml

Origin: Belarus

Raw Materials: Germany imported PVC

Packing: 500 pieces/pack





0086 13301569240



sales2@dairymachinery.cc



mobilemilkingmachine.com

Nanzhuang Village, Economic Development Zone, Yangshe Twon, Zhangjiagang City, China